

Remote Training with BrainMaster

Thomas F. Collura, Ph.D.
BrainMaster Technologies, Inc

Basics

- Session Folders
- Session Librarian
- Passkeys
- Simple Procedures

Basics

- Perform assessment
- Create plan
- Set up folders
- Convey to client
- Provide equipment to client
- Combine with passkeys to enable training

Session Folder

- Contains settings, # of sessions, duration of sessions, etc. etc.
- Contains results, summaries, statistics, etc.
- Passed back & forth between clinician and client

Session Librarian

- Use any time, any place
- “zip” up complete folder(s)
- Save on floppy disk or hard disk
- Email, ftp, SST, etc
- Folder is always complete

Passkeys

- Authorize the use of clinical sessions
- Provided to clinician by equipment provider
- Administered to client by clinician
- Ensure folders are not misused, copied, passed around, etc.

Passkeys

- Unique to device serial number
- Controlled by clinician
- Regulate session use by:
 - Time (minutes)
 - Date (any date in future)
 - Unlimited (in special cases)

Example

- Provide 40 20-minute sessions to client
- Provide 100 minutes of authorization
- Client trains 5 sessions, then system stops
- Client must report to clinician (email, phone, etc.)
- Clinician reviews results, makes changes as needed, authorizes another 100 minutes
- Etc. etc. etc.
