



BrainMaster Technologies, Inc.

From the decade of the brain into the new millenium
195 Willis Street Bedford, OH 44146 (440)232-6000 Fax: (440)232-7171
email: support@brainm.com web: www.brainm.com

MINIMUM PC COMPUTER REQUIREMENTS

BrainMaster 3.0 Software or Discovery 1.0 Software with up to 4-Channel LZT*

CPU:	Dual-Core 2.2GHz Processor
Operating System:	Windows XP, Windows Vista, Windows 7, Windows 8 Professional with Media Center Pack Add-On Windows 8.1 Professional with Media Center Pack Add-On
Memory(RAM):	Windows XP - 1GB (Minimum) 2GB (Recommended) Windows Vista - 2GB (Minimum) 4GB (Recommended) Windows 7 - 2GB (Minimum) 4GB (Recommended) Windows 8 - 2GB (Minimum) 4GB (Recommended)
Graphics Card:	512MB Dedicated OR 1GB Shared Graphics
Optical Drives:	DVD-ROM Drive: Required for BMrDVD
Input:	1 USB Port
Additional SW:	Microsoft Office: Required for Certain reports and EEGAudio Windows Media Player or 3rd Party DVD Decoder: Required for BMrDVD Adobe Acrobat Reader Adobe Flash Player: Required for BMrFlash Player Adobe Shockwave Player: Required for BMrFlash Player

Discovery 1.0 Software with BMrMMP and up to 19-Channel LZT***

CPU:	Quad-core or above
Operating System:	Windows Vista, Windows 7, Windows 8 Professional with Media Center Pack Add-On, Windows 8.1 Professional with Media Center Pack Add-On
Memory(RAM):	4GB
Graphics Card:	Direct x 10 or above compatible graphics card 1GB Dedicated
Optical Drives:	DVD-ROM Drive: Required for BMrDVD**
Input:	1 USB Port
Additional SW:	Microsoft Office: Required for certain reports and EEGAudio** Windows Media Player or 3rd Party DVD Decoder: Required for BMrDVD** Adobe Acrobat Reader Adobe Flash Player: Required for BMrFlash Player Adobe Shockwave Player: Required for BMrFlash Player

BrainAvatar™ 4.0 Software

CPU:	Quad-core or above
Operating System:	Windows Vista, Windows 7, Windows 8 Professional with Media Center Pack Add-On Windows 8.1 Professional with Media Center Pack Add-On
Memory(RAM):	4GB
Graphics Card:	Direct x 10 or above compatible graphics card 1GB Dedicated
Optical Drives:	DVD-ROM Drive: Required for BMrDVD**
Input:	1 USB Port
Additional SW:	Microsoft Office: Required for certain reports and EEGAudio** Windows Media Player or 3rd Party DVD Decoder: Required for BMrDVD** Adobe Acrobat Reader Adobe Flash Player: Required for BMrFlash Player Adobe Shockwave Player: Required for BMrFlash Player

*Live Z-Score Training is an optional purchase

**BMrMMP is an optional purchase