

## VIDEO RESOURCES FOR USE OF MINI-Q AND MINI-Q II

The instructional CD provided with the MINI-Q contains a video demonstrating the use of the original 2-channel, 6-pass MINI-Q. The relevant section of this CD begins about 2/3 of the way into the video, so you may wish to advance to that section.

This same video material is also contained in the two online Introductory Parts listed below. All of the procedures shown in this video still apply to the MINI-Q II, except for the details of the connections and selector knob. The main difference is that the MINI-Q II is able to acquire either 2 or 4 channels, and has up to 8 positions, plus a pushbutton selector for an additional 9<sup>th</sup> position.

To keep up to date on technical and applications information, you may use the BrainMaster knowledge base at:

[www.brainmaster.com/kb](http://www.brainmaster.com/kb)

And type “mini-q” into the Search: box to find current mini-q related material. This article is also online in the knowledge base, enabling you to use the links below directly.

There are online videos describing the use of the MINI-Q at:

Part I: setting up and applying electro-cap:

<http://www.youtube.com/watch?v=JBm0ryen60c>

Part II: acquiring and interpreting data:

<http://www.youtube.com/watch?v=7JaijSoMOJM>

Use of MINI-Q or MINI-Q II with z-score software:

<http://www.youtube.com/watch?v=xHASiUohrvQ>