



Minimum PC Computer Requirements

Atlantis I/II, 2EB, Discovery with BMrMMP** and up to 4 Channel LZT*

CPU:	Dual-Core 2.2GHz Processor
Operating System:	Windows XP, Windows Vista, Windows 7
Memory(RAM):	Windows XP - 1GB(Minimum) 2GB(Recommended) Windows Vista - 2GB(Minimum) 4GB(Recommended) Windows 7 - 2GB(Minimum) 4GB(Recommended)
Graphics Card:	512MB Dedicated OR 1GB Shared Graphics
Optical Drives:	DVD-ROM Drive: Required for BMrDVD**
Input:	1 USB Port
Additional SW:	Microsoft Office: Required for Certain reports and EEGAudio** Windows Media Player or 3 rd Party DVD Decoder: Required for BMrDVD** Adobe Acrobat Reader Adobe Flash Player: Required for BMrFlash Player Adobe Shockwave Player: Required for BMrFlash Player

Discovery with BMrMMP** and up to 19 Channel LZT*

CPU:	Intel® Core™ i7-740QM processor
Operating System:	Windows Vista-64 Bit, Windows 7-64 Bit
Memory(RAM):	Windows Vista-64 Bit: 6GB(Minimum) 8GB(Recommended) Windows 7-64 Bit: 6GB(Minimum) 8GB(Recommended)
Graphics Card:	1GB Dedicated
Optical Drives:	DVD-ROM Drive: Required for BMrDVD**
Input:	1 USB Port
Additional SW:	Microsoft Office: Required for certain reports and EEGAudio** Windows Media Player or 3 rd Party DVD Decoder: Required for BMrDVD** Adobe Acrobat Reader Adobe Flash Player: Required for BMrFlash Player Adobe Shockwave Player: Required for BMrFlash Player

*Live Z-Score Training is an optional purchase

**BMrMMP is an optional purchase